A Sysmex Group Company

Instructions For Use
REF: LPA 005

Prenatal 13, 18 and 21 Enumeration Probe kit

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Further information available at www.ogt.com



Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied as an essential diagnostic tool in prenatal, haematological and pathological chromosomal analysis. Target DNA, after fixation and denaturation, is available for annealing to a similarly denatured, fluorescently labelled DNA probe, which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed and the DNA is counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Information

CytoCell's prenatal Fluorescence *In Situ* Hybridisation (FISH) assays are designed for the rapid and accurate detection of the most common foetal chromosomal disorders. The probe set is intended for the detection and quantification of chromosomes 13, 18 and 21 in interphase nuclei of uncultured amniotic fluid cells by FISH. The assay should be used in conjunction with foetal karyotype analysis. Trisomy of chromosome 21, resulting in Down syndrome¹, is one of the most common chromosome abnormalities in humans and the risk of having an affected child is known to increase with maternal age². The syndrome represents a particular combination of phenotypic findings, including characteristic facial appearance, a single palmar crease and mental retardation, and may also present with hearing and heart defects. Affected individuals show a highly increased incidence of leukaemia, particularly acute megakaryocytic leukaemia³. Trisomy of chromosome 18, resulting in Edwards syndrome, occurs in around 1 in 6000-8000 live births with a female sex bias⁴. The clinical findings are variable, though many exhibit growth delay, heart defects and craniofacial anomalies, as well as possible limb and kidney abnormalities⁵. The rarest trisomy, trisomy 13, responsible for Patau syndrome occurs in approximately 1 in 16,000 newborns⁶. Individuals with Patau syndrome present abnormalities affecting many parts of the body, including the heart, spinal cord, eyes, limbs, face/skull and muscles⁶.

Probe Specification

13 unique sequence, 13q14.2 Green
18 centromere, 18p11.1 – q11.1 (D18Z1) Blue
21 unique sequence, 21q22.13 Orange

The probe set 13, 18 and 21 is a mixture of green, blue and orange directly labelled fluorescent DNA probes. The 18 centromere probe is specific for the alpha satellite DNA sequences at the D18Z1 region of the chromosome 18. The chromosome 13 probe includes RB1 gene as well as markers D13S1195, D13S1155 and D13S915. The chromosome 21 probe includes markers D21S270, D21S1867, D21S337, D21S1425, D21S1444 and D21S341.

Materials Provided

Probe: 50µl per vial (5 tests) or 100µl per vial (10 tests)
Amount of green 13q14.2 probe: 120 – 150ng/test
Amount of blue 18 centromere (D18Z1) probe: 112 – 140ng/test
Amount of orange 21q22.13 probe: 40 – 50ng/test
The probes are provided premixed in hybridisation solution (Formamide; Dextran Sulphate; SSC) and are ready to use.

Counterstain:

150µl per vial (15 tests)
The counterstain is DAPI antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide, which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The kit should be stored between -25°C to -15 °C in a freezer until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

1. Hotplate (with a solid plate and accurate temperature control up to 80°C).
2. Variable volume micropipettes and tips range 1µl - 200µl.
3. Water bath with accurate temperature control at 37°C and 72°C.
4. Microcentrifuge tubes (0.5ml).
5. Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
6. Plastic or glass coplin jars.
7. Forceps.
8. Fluorescence grade microscope lens immersion oil.
9. Bench top centrifuge.
10. Microscope slides.
11. 24x24mm coverslips.
12. Timer.
13. 37°C incubator.
14. Rubber solution glue.

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The triple bandpass filter DAPI/FITC/TRITC is optimal for viewing the green and orange fluorophores as well as counterstain simultaneously. The blue fluorophore has specificity to the Aqua and DEAC spectrum (single bandpass Aqua or DEAC filter is required). The triple bandpass filter DAPI/FITC/Texas Red can be used to view all three fluorophores and DAPI simultaneously.

Sample Preparation

The prenatal kit is designed for use on uncultured amniocytes fixed in Carnoy's fixative (see procedure below). Amniotic fluid sample collection should be performed according to the laboratory or institution guidelines. Amniotic fluid samples that appear bloody or brown should not be used, since they may contain maternal blood and may lead to false results.

Suggested Protocol

Preparation of fresh amniotic fluid samples for FISH:

1. Centrifuge 2 – 5ml of whole amniotic fluid specimen for 7 minutes at 180xg, carefully remove the supernatant without disturbing the cell pellet.
2. Resuspend the pellet in 2ml of 0.075M KCl. Leave at room temperature (RT) for 5 minutes.
3. Add 2ml of fresh fixative (3:1 methanol:glacial acetic acid) to the cells/hypotonic solution, adding the first ml dropwise whilst continuously mixing. Mix well.
4. Centrifuge the suspension for 5 minutes at 280xg, carefully remove the supernatant and resuspend the pellet in 2ml of fresh fixative.
5. Fixed specimens can be stored at this stage in a freezer at -20°C.
6. If the sample is not to be frozen, centrifuge the tube at 280xg for 5 minutes. Remove as much supernatant as possible without disturbing the cell pellet. Flick the tube to resuspend the pellet in the small amount of fluid remaining.
7. To prepare slides for FISH, spot the cell suspension directly onto slide. Allow to air dry.

Recommended slide Pretreatment:

1. Immerse the slide prepared from uncultured amniocytes in 2xSSC for 1 hour at 37°C.
2. Place the slide in freshly made pepsin working solution (5mg of pepsin added to 100ml of 0.01M HCl) for 13 minutes at 37°C.
3. Immerse the slide in phosphate buffered saline (PBS) at RT for 5 minutes.
4. Immerse the slide in post fixation solution (0.95% formaldehyde: 1.0ml of 37% formaldehyde, 0.18g of MgCl₂ and 39.0ml of PBS) for 5 minutes at RT.
5. Immerse the slide in PBS at RT for 5 minutes.
6. Immerse the slide in 70% ethanol at RT. Allow the slide to stand in the ethanol wash for 1 minute.
7. Remove the slide from 70% ethanol. Repeat step 6 with 85% ethanol, followed by 100% ethanol.
8. Allow to air dry.

FISH Protocol

(Note: Please ensure that exposure of the probe to laboratory lights is limited at all times)

Slide preparation (skip this step if the slide was pretreated according to the protocol above)

1. Spot the cell sample onto a glass microscope slide. Allow to dry.
2. Immerse the slide in 2xSSC for 2 minutes at RT without agitation.
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 minutes at RT.
4. Allow to dry.

Pre-Denaturation

5. Remove the probe from the freezer and allow it to warm to RT.
6. Ensure that the probe solution is uniformly mixed with a pipette.
7. Remove 10µl of probe per test, and transfer it to a microcentrifuge tube. Quickly return the remaining probe to the freezer.
8. Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.
9. Spot 10µl of probe mixture onto the cell sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

Denaturation

10. Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 2 minutes.

Hybridisation

11. Place the slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

Post-Hybridisation Washes

12. Remove the coverslip and all traces of glue carefully.
13. Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
14. Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.
15. Drain the slide and apply 10µl of DAPI antifade onto each hybridisation area.
16. Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
17. View with a fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below RT.

Procedural Recommendations

1. Baking or ageing of slides is not recommended as it may reduce signal fluorescence.
2. Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by CytoCell Ltd.
3. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators as these temperatures are critical for optimum product performance.
4. The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
5. Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.

Interpretation of Results

The sensitivity and specificity of FISH depends on a number of parameters, which vary from one cell type to another; from one probe to another; with the cell techniques used; and within the individual laboratory. Therefore we recommend for the use of the prenatal kit, that each laboratory should have its own standard material and determine its own FISH assay cutoff values for karyotypically normal and aneuploid samples (for guidelines contact CytoCell).

Expected Results

A normal cell should show 2 green, 2 blue and 2 orange signals (2G, 2B, 2O). Cells with trisomy 13 should show 3 green, 2 blue and 2 orange signals (3G, 2B, 2O). Cells with trisomy 18 should show 2 green, 3 blue and 2 orange signals (2G, 3B, 2O). Cells with trisomy 21 should show 2 green, 2 blue and 3 orange signals (2G, 2B, 3O).

Limitations

This test does not detect structural chromosome abnormalities and mosaicism. It also does not detect numerical abnormalities of chromosomes other than tested. The test is not designed to evaluate the risk of trisomy. Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information.

This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone.

Additional Information

For additional product information please contact the CytoCell Technical Support Department.

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Patents and Trademarks
CytoCell is a registered trademark of CytoCell Ltd.