A Sysmex Group Company

Instructions For Use  
REF: LPS 047-S / LPS 047

## 1p36/1q25 and 19q13/19p13 Deletion Probe Kit

PROFESSIONAL USE ONLY

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Further information available at [www.ogt.com](http://www.ogt.com)

### Intended Use

The 1p36/1q25 and 19q13/19p13 Deletion Probe Kit is intended to identify, via Fluorescence *In Situ* Hybridisation (FISH), deletions of the 1p36.32 region on chromosome 1 and deletions of the 19q13.33 region on chromosome 19 in patients with glial brain tumours. It is intended to be used on formalin fixed paraffin embedded (FFPE) tissue samples. This product is for professional use only and is intended to be an adjunct to clinical cytogenetics.

### Background

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied to the assessment of solid tumour biopsies as well, which can provide important information for prediction of tumour progression. Current methodologies, namely immunohistochemistry or blotting can provide information at the gene expression level but, when carried out on tissue sections (either cryostat or paraffin embedded), FISH can provide information at a gene level, *in situ*, at the precise site within the tumour. This can reveal cell-to-cell heterogeneity and enable the detection of small clones of genetically distinct cells.

### Probe Information


Deletions of the 1p36.32 region including the TP73 (*tumor protein 73*) gene and deletions of the 19q13.33 region including the GLTSCR1 and GLTSCR2 (*glioma tumor suppressor candidate region genes 1 and 2*) genes are frequently reported in cases of glial tumours.

Astrocytomas and oligodendrogliomas are the most common gliomas that arise from glial cells. They make up about 40% of all CNS tumours<sup>1</sup> and more than 60% of primary brain cancers<sup>2</sup>.

Concurrent losses, 'co-deletion', of the 1p36.32 and 19q13.33 regions are reported in approximately 80% of oligodendrogliomas, two-thirds of anaplastic oligodendrogliomas, as well as subsets of oligoastrocytomas and anaplastic oligoastrocytomas<sup>3,4</sup>; the majority of these losses have been shown to be mediated by the presence of an unbalanced t(1;19)(q10;p10) translocation. The presence of a 1p and 19q co-deletion is a strong prognostic factor in these diseases, where it is associated with improved prognosis and responsiveness to therapy<sup>5</sup>. 1p and 19q co-deletion has also been shown to occur in a subset of extraventricular neurocytomas, and may be associated with aggressive histology in these tumours<sup>6</sup>.

### Probe Specification


**Probe 1:**  
1p36.32, Red  
1q25.2 Green

The 1p36.32 probe, labelled in red, is 433kb in size and covers the region between markers RH122382 and SHGC-74088. The 1q25.2 probe, labelled in green,

consists of three probes (143kb, 334kb and 140kb) that cover regions including markers SHGC-75984 and SHGC-147545.

**Probe 2:**  
19p13.2 Green  
19q13.33, Red

The 19p13.2 probe, labelled in green, consists of three probes (148kb, 174kb and 131kb) that cover regions including markers D19S1025 and D19S677E. The 19q13.33 probe, labelled in red, is 357kb in size and covers the region between markers RH91552 and D19S902.

### Materials Provided

50µl per vial (5 tests) or 100µl per vial (10 tests)

#### Probe 1:

Amount of red 1p36.32 probe: 80-100ng/test  
Amount of green 1q25.2 probe: 240-300ng/test

#### Probe 2:

Amount of red 19q13.33 probe: 60-75ng/test  
Amount of green 19p13.2 probe: 264-330ng/test

The probes are provided premixed in hybridisation solution (Formamide; Dextran Sulphate; SSC) and are ready to use.

**Counterstain:** 150µl per vial (15 tests) or 300µl per vial (30 tests)

The counterstain is DAPI antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)).

### Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide, which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. Dispose of all hazardous materials according to your institution's guidelines for hazardous waste disposal.
6. Operators must be capable of visually distinguishing between red, blue and green.

### Storage and Handling

The kit should be stored between -25°C to -15 °C in a freezer until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

### Equipment Necessary but not Supplied

1. Hotplate (with a solid plate and accurate temperature control up to 80°C).
2. Variable volume micropipettes and tips range 1µl - 200µl.
3. Water bath with accurate temperature control at 72°C.
4. Microcentrifuge tubes (0.5ml).
5. Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
6. Plastic or glass coplin jars.
7. Forceps.
8. Fluorescence grade microscope lens immersion oil.
9. Bench-top centrifuge.
10. Microscope slides.
11. 24x24mm coverslips.
12. Timer.
13. 37°C incubator.
14. Rubber solution glue.
15. Tissue Pretreatment Kit (LPS 100).

### Fluorescence Microscope Recommendation

Use a 100-watt mercury lamp and plan apochromat objectives x63 or x100 for optimal visualisation. Use a triple bandpass filter DAPI/FITC/Texas Red for optimal visualisation of the green and red fluorophores and DAPI simultaneously. Check the fluorescence microscope before use to ensure it is operating correctly. Use immersion oil that is suitable for fluorescence microscopy and formulated for low autofluorescence. Follow manufacturers' recommendations in regards to the life of the lamp and the age of the filters.

### Sample Preparation

The kit is designed for use on formalin fixed paraffin embedded (FFPE) tissue sections or tissue microarrays (TMA) and should be prepared according to the laboratory or institution guidelines. For FISH, 4µm - 6µm thick FFPE tissue sections should be used.

## Tissue Sample Pretreatment

Tissue sample pretreatment should be done according to the laboratory or institution guideline. For optimal results use the Tissue Pretreatment Kit (LPS 100).

## FISH Protocol

(Note: Ensure that exposure of the probe to laboratory lights is limited at all times)

### Pre-Denaturation

1. Remove the probe from the freezer and allow it to warm to room temperature (RT). Briefly centrifuge tubes before use
2. Ensure that the probe solution is uniformly mixed with a pipette.
3. Remove 10µl -15µl of probe per test, and transfer it to a microcentrifuge tube. Quickly return the remaining probe to the freezer.
4. Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.
5. Spot 10µl - 15µl of probe mixture onto the sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

### Denaturation

6. Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 5 minutes.

### Hybridisation

7. Place the slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

### Post-Hybridisation Washes

8. Remove the coverslip and all traces of glue carefully.
9. Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
10. Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.
11. Drain the slide and apply 10µl - 15µl of DAPI antifade onto each sample.
12. Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
13. View with a fluorescence microscope.

### Comments

Hybridisation efficiency and tissue morphology are usually negatively correlated. Aggressive pretreatment procedures improving hybridisation efficiency (e.g. an extended enzyme digestion time) tend to destroy cell structure and tissue morphology. However, mild pretreatment saving tissue structures may not be sufficient for probe penetration and successful FISH results.

The optimal length of heat pretreatment and enzyme digestion time will depend on the age of the block, the tissue composition, and quality of tissue fixation. Enzyme digestion should be decreased for core biopsies and any sections that contain few tumour cells or have large areas of necrosis. These samples need to be handled particularly carefully to avoid over-digestion.

### Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark below 4°C.

### Procedural Recommendations

1. Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by CytoCell Ltd.
2. Use a calibrated thermometer for measuring temperatures of solutions, waterbaths and incubators as these temperatures are critical for optimum product performance.
3. The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
4. Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.
5. Over hybridisation can result in additional or unexpected signals.
6. Users should optimise the protocol for their own samples prior to using the test for diagnostic purposes.

### Expected Results

#### Probe 1:

A normal cell should show 2 red and 2 green signals (2R, 2G). Cells carrying a 1p36 deletion should show 1 red and 2 green signals (1R, 2G).

#### Probe 2:

A normal cell should show 2 red and 2 green signals (2R, 2G). Cells carrying a 19q13 deletion should show 1 red and 2 green signals (1R, 2G).

### Known Cross-Reactivity

No known cross-reactivity.

### Limitations

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone.

Failure to adhere to the protocol may alter the performance of the assay and may yield erroneous results.

This kit has not been validated for purposes outside of the intended use stated.

This device is designed to detect genomic losses larger than the regions covered by the red clones in this probe set, which include the 1p36.32 and 19q13.33 regions. Genomic losses outside these regions or partial losses of these regions may not be detected with this product.

### Additional Information

For additional product information please contact the CytoCell Technical Support Department.

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

W: www.ogt.com

### References/Bibliographie/Bibliographia/Literatur/Bibliografía

1. CBTRUS (2004). Statistical Report: Primary Brain Tumors in the United States, 1997-2001. Central Brain Tumor Registry of the United States.
2. Thompson L. World Health Organization classification of tumours: pathology and genetics of head and neck tumours. Ear Nose Throat J. 2006 Feb;85(2):74.
3. Vogazianou AP, Chan R, Bäcklund LM, Pearson DM, Liu L, Langford CF, et al. Distinct patterns of 1p and 19q alterations identify subtypes of human gliomas that have different prognoses. Neuro Oncol. 2010 Jul;12(7):664-78.
4. Bromberg JEC, et al., Oncol. 2009. 14:155-163.
5. Jenkins RB et al., Cancer Res 2006;66(20):9852-61.
6. Rodriguez FJ et al., Brain Pathol 2009;19(4):623-9.

	<b>EN:</b> Catalogue number <b>DE:</b> Bestellnummer <b>FR:</b> Référence du catalogue <b>IT:</b> Numero di catalogo <b>ES:</b> Número de catálogo
	<b>EN:</b> <i>In vitro</i> diagnostic device <b>DE:</b> <i>In-vitro</i> -Diagnostikum <b>FR:</b> Dispositif médical de diagnostic <i>in vitro</i> <b>IT:</b> Dispositivo medico-diagnostico <i>in vitro</i> <b>ES:</b> Producto sanitario para diagnóstico <i>in vitro</i>
	<b>EN:</b> Batch code <b>DE:</b> Loscode <b>FR:</b> Code du lot <b>IT:</b> Codice lotto <b>ES:</b> Código de lote
	<b>EN:</b> Consult instructions for use <b>DE:</b> Gebrauchsanweisung beachten <b>FR:</b> Consulter la notice d'utilisation <b>IT:</b> Consultare le istruzioni per l'uso <b>ES:</b> Consulte las instrucciones de uso
	<b>EN:</b> Manufacturer <b>DE:</b> Hersteller <b>FR:</b> Fabricant <b>IT:</b> Prodotto da <b>ES:</b> Fabricante
	<b>EN:</b> Use by <b>DE:</b> Haltbarkeitsdatum <b>FR:</b> A utiliser avant <b>IT:</b> Scadenza <b>ES:</b> Fecha de caducidad
	<b>EN:</b> Temperature limitation <b>DE:</b> Temperaturbegrenzung <b>FR:</b> Limites de température <b>IT:</b> Limiti di temperatura <b>ES:</b> Limitación de temperatura
	<b>EN:</b> Sufficient for <n> tests <b>DE:</b> Ausreichend für <n> Tests <b>FR:</b> Suffisant pour <n> tests <b>IT:</b> Sufficiente per <n> test <b>ES:</b> Válido para <n> análisis
	<b>EN:</b> Contents <b>DE:</b> Inhalt <b>FR:</b> Contenu <b>IT:</b> Contenuto <b>ES:</b> Contenido

### Patents and Trademarks

CytoCell is a registered trademark of Cytozell Ltd.

This product contains technology licensed from Life Technologies Corporation and is available for human diagnostics or life science research use only.

#### Cytocell Ltd.

Oxford Gene Technology,  
418 Cambridge Science Park,  
Milton Road,  
Cambridge, CB4 0PZ, UK  
T: +44(0)1223 294048  
F: +44(0)1223 294986  
E: probes@cytozell.com  
W: www.ogt.com