A Sysmex Group Company

Instructions For Use

REF: LPE 008B / LPE 012B / LPE 017B

Satellite Enumeration Probes

PROFESSIONAL USE ONLY

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Further information available at www.ogt.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied as an essential diagnostic tool in prenatal, haematological and pathological chromosomal analysis. Target DNA, after fixation and denaturation, is available for annealing to a similarly denatured, fluorescently labelled DNA probe, which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed and the DNA is counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Limitations

The test is not intended for: use as a stand-alone diagnostic, prenatal testing, population-based screening, near-patient testing or self-testing. This product is intended for laboratory professional use only; all results should be interpreted by suitably-qualified staff, taking into account other relevant test results.

This product has not been validated for use on sample types or disease types other than those specified in the intended use.

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone.

Failure to adhere to the protocol may affect the performance and lead to false positive/negative results.

This kit has not been validated for purposes outside of the intended use stated.

Probe Information

Satellite probes are specific for human chromosomes. They are highly repeated human DNA sequences found in the centromere, pericentromeric or heterochromatic block of each of the 24 chromosomes. The probes enable the identification and enumeration of human chromosomes in interphase cells or metaphase chromosomes from peripheral blood samples.

Probe Specification

The probes are produced in a concentrated form to allow mixing, if required, of up to three probes in the same hybridisation, from CytoCell's range of concentrated Satellite probes. A final volume of 10µl of probe solution is required per hybridisation.

The probes are directly labelled with a blue fluorophore (Aqua or DEAC spectrum). For detailed probe specifications refer to Table 1.

Table 1: Probe Specifications

Chr	Catalogue Number*	Locus	Chromosome Region	DNA Class
8	LPE 008B	D8Z2	8p11.1-q11.1	α-satellite
12	LPE 012B	D12Z3	12p11.1-q11.1	α-satellite
17	LPE 017B	D17Z1	17p11.1-q11.1	α-satellite

*B specifies a blue label

This kit contains only one of the probes from the range of directly labelled blue human alpha satellite probes.

Materials Provided

Probe: 30µl per vial (10 tests)

Amount of blue D8Z2 probe: 63.5-91.5ng/test

Amount of blue D12Z3 probe: 19.2ng/test

Amount of blue D17Z1 probe: 14.4ng/test

The probe is produced in a concentrated form. It is provided in hybridisation solution (Formamide; Dextran Sulphate; SSC).

Hybridisation solution (Formamide; Dextran Sulphate; SSC): 150µl per vial

Counterstain: 150µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

- For *in vitro* diagnostic use. For professional use only.
- Wear gloves when handling DNA probes and DAPI counterstain.
- Probe mixtures contain formamide, which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
- DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
- Dispose of all hazardous materials according to your institution's guidelines for hazardous waste disposal.
- Operators must be capable of visually distinguishing between red, blue and green.

Storage and Handling

The kit should be stored between -25°C to -15 °C in a freezer until the expiry date indicated on the kit label. Store the probe and counterstain vials in the dark. Ensure that exposure of the probe and counterstain to laboratory lights is limited at all times.

Equipment Necessary but not Supplied

- Hotplate (with a solid plate and accurate temperature control up to 80°C).
- Variable volume micropipettes and tips range 1µl - 200µl.
- Water bath with accurate temperature control at 72°C.
- Microcentrifuge tubes (0.5ml).
- Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
- Plastic or glass coplin jars.
- Forceps.
- Fluorescence grade microscope lens immersion oil.
- Bench top centrifuge.
- Microscope slides.
- 24x24mm coverslips.
- Timer.
- 37°C incubator.
- Rubber solution glue.

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100-watt mercury lamp and plan apochromat objectives x63 or x100. The blue fluorophore has specificity to the Aqua and DEAC spectrum (single bandpass Aqua or DEAC filter is required).

Check the fluorescence microscope before use to ensure it is operating correctly. Use immersion oil that is suitable for fluorescence microscopy and formulated for low autofluorescence. Follow manufacturers' recommendations in regards to the life of the lamp and the age of the filters.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative, that should be prepared according to the laboratory or institution guidelines. Prepare air dried samples on microscope slides according to standard cytogenetic procedures.

FISH Protocol

(Note: Ensure that exposure of the probe and counterstain to laboratory lights is limited at all times).

Slide preparation

- Spot the cell sample onto a glass microscope slide. Allow to dry.
- Immerse the slide in 2xSSC for 2 minutes at room temperature (RT) without agitation.
- Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 minutes at RT.
- Allow to dry.

Pre-Denaturation

- Remove the probe from the freezer and allow it to warm to RT. Briefly centrifuge tubes before use.
- Ensure that the probe solution is uniformly mixed with a pipette.
- Using fresh pipette tips remove (final volume of 10µl of probe solution):
 - for a **single probe hybridisation**: 3µl of probe and 7µl of hybridisation solution per test
 - for a **two probe hybridisation**: 3µl of each probe and 4µl of hybridisation solution per test
 - for a **three probe hybridisation**: 3µl of each probe and 1µl of hybridisation solution per test
- and transfer it to a microcentrifuge tube, gently vortex to mix and pulse spin in a microcentrifuge. Quickly return the remaining probe to the freezer.
- Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.

9. Spot 10µl of probe mixture onto the cell sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

Denaturation

10. Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 2 minutes.

Hybridisation

11. Place the slide in a humid, lightproof container at 37°C (+/- 1°C) for 1 hour to overnight.

Post-Hybridisation Washes

12. Remove the coverslip and all traces of glue carefully.
13. Immerse the slide in 0.25xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation*.
14. Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.
15. Drain the slide and apply 10µl of DAPI antifade onto each sample.
16. Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
17. View with a fluorescence microscope.

* If final signal is poor, repeat FISH using 0.4xSSC post-hybridisation wash.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at/or below RT.

Procedural Recommendations

1. Baking or ageing of slides may reduce signal fluorescence.
2. Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by CytoCell Ltd.
3. Use a calibrated thermometer for measuring temperatures of solutions, waterbaths and incubators as these temperatures are critical for optimum product performance.
4. The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
5. Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.
6. Over hybridisation can result in additional or unexpected signals.
7. Users should optimise the protocol for their own samples prior to using the test for diagnostic purposes.

Expected Results for a single probe hybridisation

There will be differences in the relative size of signals observed between chromosomes due to the difference in copy number of repeat sequences between chromosomes. A diploid sample should show a fluorescent signal at the centromere of both of the corresponding chromosomes in 70 – 90% of cells analysed.

Known Cross-Reactivity

The LPE017B probe may show faint cross hybridisation to the centromeric region of chromosomes 11 and X.

The LPE012B may show faint cross-hybridisation to the centromeric region of chromosomes 3, 6, 7 and 10.

Limitations

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone.
Failure to adhere to the protocol may affect the performance and results.

Additional Information

For additional product information please contact the CytoCell Technical Support Department.
T: +44 (0)1223 294048
E: techsupport@cytozell.com
W: www.cgt.com

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Patents and Trademarks
CytoCell is a registered trademark of CytoCell Ltd.