A Sysmex Group Company

Instructions For Use
REF: LPS 046-S / LPS 046

ROS1 Plus Breakapart Probe

PROFESSIONAL USE ONLY

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Further information available at www.ogt.com

Intended Use

The ROS1 Plus Breakapart Probe is intended to identify, via fluorescence *in situ* hybridisation (FISH), rearrangements or deletions of the ROS1 region on chromosome 6 at location q22.1. It is intended to be used on formalin-fixed paraffin-embedded (FFPE) tissues or tissue microarrays (TMAs). This product is for professional use only and is intended to be an adjunct to Clinical Cytogenetics.

Background

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied to the assessment of solid tumour biopsies as well, which can provide important information for prediction of tumour progression. Current methodologies, namely immunohistochemistry or blotting can provide information at the gene expression level but, when carried out on tissue sections (either cryostat or paraffin embedded), FISH can provide information at a gene level, *in situ*, at the precise site within the tumour. This can reveal cell-to-cell heterogeneity and enable the detection of small clones of genetically distinct cells.


Probe Information

The ROS1 (*ROS proto-oncogene 1, receptor tyrosine kinase*) gene at 6q22 is an ALK (anaplastic lymphoma receptor tyrosine kinase) gene paralogue which encodes a type I integral membrane protein with tyrosine kinase activity¹. ROS1 rearrangements define a molecular subset of non-small cell lung cancer (NSCLC) and are seen in approximately 2% of patients with NSCLC². A number of partner genes have been identified, including SLC34A2, CD74 and SDC4³. It has been shown that these ROS1 fusions activate the pSTAT3, PI3K/AKT/mTOR and SHP-2 phosphatase pathways^{4,5}. NSCLC patients with ROS1 rearrangements have been shown to respond to treatment with ALK/MET tyrosine kinase inhibitors, such as crizotinib⁶. ROS1 rearrangements with the GOPC (*golgi associated PDZ and coiled-coil motif containing*) gene fusion partner were originally reported in glioblastoma, but have now also been detected on cholangiocarcinoma, ovarian cancer and NSCLC patient samples⁷⁻⁹.

The ROS1 Plus design covers the ROS1 region and the region deleted in ROS1-GOPC fusions.

Probe Specification

ROS1 Plus, 6q22.1, Green
ROS1 Plus, 6q22.1, Red

The ROS1 Plus probe consists of 3 probes (175kb, 123kb, 184kb), labelled in green, situated proximal to the ROS1 gene and including markers D6S2428 and

SHGC-149552 and a red probe (206kb) covering the distal part of the ROS1 gene and a region up to RH78348.

Materials Provided

Probe: 50µl per vial (5 tests) or 100µl per vial (10 tests)

Amount of green ROS1 Plus probe: 411-615ng/test

Amount of red ROS1 Plus probe: 102-177ng/test

The probes are provided premixed in hybridisation solution (Formamide; Dextran Sulphate; SSC) and are ready to use.

Counterstain: 150µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide, which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. Dispose of all hazardous materials according to your institution's guidelines for hazardous waste disposal.
6. Operators must be capable of visually distinguishing between red, blue and green.

Storage and Handling

The kit should be stored between -25°C to -15°C in a freezer until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

1. Hotplate (with a solid plate and accurate temperature control up to 80°C).
2. Variable volume micropipettes and tips range 1µl - 200µl.
3. Water bath with accurate temperature control at 72°C.
4. Microcentrifuge tubes (0.5ml).
5. Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
6. Plastic or glass coplin jars.
7. Forceps.
8. Fluorescence grade microscope lens immersion oil.
9. Bench-top centrifuge.
10. Microscope slides.
11. 24x24mm coverslips.
12. Timer.
13. 37°C incubator.
14. Rubber solution glue.
15. Tissue Pretreatment Kit (LPS 100).

Fluorescence Microscope Recommendation

Use a 100-watt mercury lamp and plan apochromat objectives x63 or x100 for optimal visualisation. Use a triple bandpass filter DAPI/FITC/Texas Red for optimal visualisation of the green and red fluorophores and DAPI simultaneously. Alternatively for viewing red and green fluorophores use dual bandpass filter FITC/Texas Red.

Check the fluorescence microscope before use to ensure it is operating correctly. Use immersion oil that is suitable for fluorescence microscopy and formulated for low autofluorescence. Follow manufacturers' recommendations in regards to the life of the lamp and the age of the filters.

Sample Preparation

The kit is designed for use on Formalin Fixed Paraffin Embedded (FFPE) tissue sections or Tissue Microarrays (TMA) and should be prepared according to the laboratory or institution guidelines. Use 3µm - 5µm thick FFPE tissue sections for FISH.

Tissue Sample Pretreatment

Tissue sample pretreatment should be done according to the laboratory or institution guideline. Use the Tissue Pretreatment Kit (LPS 100) for optimal results.

FISH Protocol

(Note: Ensure that exposure of the probe to laboratory lights is limited at all times)

Pre-Denaturation

1. Remove the probe from the freezer and allow it to warm to room temperature (RT). Briefly centrifuge tubes before use
2. Ensure that the probe solution is uniformly mixed with a pipette.
3. Remove 10µl - 15µl of probe per test, and transfer it to a microcentrifuge tube. Quickly return the remaining probe to the freezer.
4. Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.
5. Spot 10µl - 15µl of probe mixture onto the sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

Denaturation

6. Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 5 minutes.

Hybridisation

7. Place the slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

Post-Hybridisation Washes

8. Remove the coverslip and all traces of glue carefully.
9. Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
10. Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.

DS280/CE v006.00/2020-12-01 (S060 v3)

Page 1 of 6

- Drain the slide and apply 10µl - 15µl of DAPI antifade onto each sample.
- Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
- View with a fluorescence microscope.

Comments

Hybridisation efficiency and tissue morphology are usually negatively correlated. Aggressive pretreatment procedures improving hybridisation efficiency (e.g. an extended enzyme digestion time) tend to destroy cell structure and tissue morphology. However, mild pretreatment saving tissue structures may not be sufficient for probe penetration and successful FISH results. The optimal length of heat pretreatment and enzyme digestion time will depend on the age of the block, the tissue composition, and quality of tissue fixation. Decrease enzyme digestion for core biopsies and any sections that contain few tumour cells or have large areas of necrosis. Handle these samples particularly carefully to avoid over-digestion.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark below 4°C.

Procedural Recommendations

- Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by CytoCell Ltd.
- Use a calibrated thermometer for measuring temperatures of solutions, waterbaths and incubators as these temperatures are critical for optimum product performance.
- The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
- Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.
- Over hybridisation can result in additional or unexpected signals.
- Users should optimise the protocol for their own samples prior to using the test for diagnostic purposes.

Expected Results

In a normal cell, two red/green fusion signals are expected (2F). In a cell with a balanced ROS1 rearrangement, the expected signal pattern will be 1 red, 1 green and 1 fusion (1R, 1G, 1F). In the event of a deletion of material distal to ROS1 there will be loss of one red signal, and the expected signal pattern will be 1 green and 1 fusion (1G, 1F). Other signal patterns are possible in aneuploid/unbalanced specimens.

Note: The proximal and distal probes in LPS 046, ROS1 Plus Breakapart probe, are labelled in the opposite colours to the proximal and distal probes in CytoCell LPS 022, ROS1 Breakapart probe.

Known Cross-Reactivity

No known cross-reactivity.

Limitations

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone.

Failure to adhere to the protocol may affect the performance and results. This kit has not been validated for purposes outside of the intended use stated.

Additional Information

For additional product information please contact the CytoCell Technical Support Department.

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Patents and Trademarks
CytoCell is a registered trademark of CytoCell Ltd.